



germanengineering

In Deutschland: Vertrieb ausschließlich
über medizinischen Fachhandel

BC 30 Blutdruckmessgerät

Vollautomatische Blutdruck- und Pulsmessung
am Handgelenk

Übersichtliches LCD-Display

Durchschnittswert der letzten 3 gespeicherten Messwerte

2 x 60 Speicherplätze

Mit WHO*-Indikator (Einstufung der Messergebnisse gemäß
den WHO*-Richtlinien)

Arrhythmie-Erkennung: erkennt und warnt bei eventuellen
Herzrhythmusstörungen

Medizinprodukt

Manschettengröße für Handgelenksumfänge
von 13,5 - 19,5 cm

Datum und Uhrzeit /Abschaltautomatik

Meldung bei Anwendungsfehler

Automatische Druckvorwahl und Luftablass

Blutdruckpass

Batteriewechselanzeige

Praktische Aufbewahrungsbox

Inkl. 2 x 1,5V Batterien AAA

Produktmaße (L/B/H): 72 x 70 x 27 mm

Produktgewicht: ca. 111 g

3 Jahre Garantie

VE: 4 / Transportkarton: 8 x 4

EAN-Nr.: 4211125 65051 3

Art.-Nr.: 650.51

BC 30 Blood pressure monitor

Fully automatic blood pressure and pulse measurement
on wrist

Clearly arranged LCD display

Average of the last 3 saved measured values

2 x 60 memory spaces

With WHO* indicator (classification of values according
to WHO* guidelines)

Arrhythmia detection: warning in case of possible
heart rhythm disturbance

Medical product

Cuff size for wrist circumference 13.5 - 19.5 cm

Date and time / Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Blood pressure passport

Low battery indicator

Practical storage box

Incl. 2 x 1.5V AAA batteries

Product measurements (LW/H): 72 x 70 x 27 mm

Product weight: approx. 111 g

3 year guarantee

Sales unit: 4 / Shipping carton: 8 x 4

EAN no.: 4211125 65051 3

Item no.: 650.51

*WHO = World Health Organization (Weltgesundheitsorganisation)